

ИНФУЗОМАТ

ШПРИЦ 50 МЛ	СКОРОСТЬ МЛ/ЧАС							
ДОФАМИН 200 МГ МКГ/КГ/МИН	КГ 50	КГ 60	КГ 70	КГ 80	КГ 90	КГ 100	КГ 110	КГ 120
1	0.8	0.9	1.1	1.2	1.4	1.5	1.7	1.8
2	1.5	1.8	2.1	2.4	2.7	3.0	3.3	3.6
3	2.3	2.7	3.2	3.6	4.1	4.5	5.0	5.4
4	3.0	3.6	4.2	4.8	5.4	6.0	6.6	7.2
5	3.8	4.5	5.3	6.0	6.8	7.5	8.3	9.0
6	4.5	5.4	6.3	7.2	8.1	9.0	9.9	10.1
7	5.3	6.3	7.4	8.4	9.5	10.5	11.6	12.6
8	6.0	7.2	8.4	9.6	10.8	12.0	13.2	14.4
9	6.8	8.1	9.5	10.8	12.2	13.5	14.9	16.2
10	7.5	9.0	10.5	12.0	13.5	15.0	16.5	18.0
11	8.3	9.9	11.6	13.2	14.9	16.5	18.2	19.8
12	9.0	10.8	12.6	14.4	16.2	18.0	19.8	21.6
13	9.8	11.7	13.7	15.6	17.6	19.5	21.5	23.4
14	10.5	12.6	14.7	16.8	18.9	21.0	23.1	25.2
15	11.3	13.5	15.8	18.0	20.3	22.5	24.8	27.0
16	12.0	14.4	16.8	19.2	21.6	24.0	26.4	28.8
17	12.8	15.3	17.9	20.4	23.0	25.5	28.1	30.6
18	13.5	16.2	18.9	21.6	24.3	27.0	29.7	32.4
19	14.3	17.1	20.0	22.8	25.7	28.5	31.4	34.2
20	15.0	18.0	21.0	24.0	27.0	30.0	33.0	36.0